



# Blood Collection Tubes

Vacuum - VAC-T / Non Vacuum - RAP-T



**VAC-T**



**RAP-T**



Mktd. By:  
**RAPID DIAGNOSTIC PVT. LTD.**  
ISO 9001 : 2000 CERTIFIED COMPANY



Local Contact:



**VAC-T / RAP-T**

Tube Type & Color Code	Tube (mm)	Additive	No. of Inversions (Mixing)	Draw Volume	Laboratory Use
------------------------	-----------	----------	----------------------------	-------------	----------------

**K<sub>2</sub>/K<sub>3</sub> EDTA**

 Lavender (Purple)	13x75	K <sub>2</sub> EDTA	5-10	2.0ml	Whole Blood, Hematology, Blood Banking
 Lavender (Purple)	13x75	K <sub>3</sub> EDTA	5-10	2/3.0ml	Whole Blood, Hematology, Blood Banking


**PLAIN TUBE WITH CLOT ACTIVATOR**

 Red	13x75	Clot Activator	2-3	5.0ml	Biochemistry, Serology, Blood Banking
--	-------	----------------	-----	-------	---------------------------------------

**SERUM SEPARATING TUBE (SST) WITH GEL SEPARATOR**

 Golden Yellow	13x75	Gel Separator	5-8	4.0ml	Biochemistry, Serology, Blood Banking
--	-------	---------------	-----	-------	---------------------------------------

**SODIUM FLUORIDE**

 Grey	13x75	Sodium Fluoride	5-10	2.0ml	Glucose/Alcohol
---	-------	-----------------	------	-------	-----------------

**SODIUM CITRATE TUBE**

 Blue	13x75	3.2% Sodium Citrate	5-8	1.8ml	Coagulation/ Factor Assay
---	-------	---------------------	-----	-------	---------------------------

**VAC-T / RAP-T**

Tube Type & Color Code	Tube (mm)	Additive	No. of Inversions (Mixing)	Draw Volume	Laboratory Use
------------------------	-----------	----------	----------------------------	-------------	----------------

**SODIUM HEPARIN TUBE**

 Green	13x75	Sodium Heparin	5-10	5.0ml	Hematology/ Chemistry (Plasma) Assay
--	-------	----------------	------	-------	--------------------------------------

**ESR TUBE**

 Black	13x75	3.8% Sodium Citrate	2-3	1.6ml	ESR
--	-------	---------------------	-----	-------	-----

**NEEDLE**



**NEEDLE HOLDER**



**ESR TUBE FOR ESR ANALYZERS**

**Glass tube**



ESR - T	8x120mm	3.8% Sodium Citrate	2-3	1.6ml	ESR
---------	---------	---------------------	-----	-------	-----